

## HPLC CHROMATOGRAPY – TOPIC TEST 1

### QUESTION 1

In HPLC systems:

- A The mobile phase can be a liquid or a gas.
- B The stationary phase is inside a metal column.
- C The sample needs to be vaporised before it can be injected.
- D A sample cannot be recovered as it is destroyed.

### QUESTION 2

In HPLC, a substance with a high retention time:

- A Has a high degree of adsorption into the stationary phase.
- B Has a high degree of adsorption into the mobile phase.
- C Remains in solution for long periods of time.
- D Takes a long time to go through the detector.


### QUESTION 3

Which of the following is not true about HPLC systems:

- A The mobile phase is pumped under pressure.
- B A variety of column packing materials is available.
- C It is suitable for both heat resistant and heat sensitive substances.
- D It is widely used in food analysis.

### QUESTION 4


The following represents a schematic HPLC system:



- A The sample is injected at Z.
- B The solvent is stored at Y.
- C The sample is recovered at Y.
- D The chromatogram is generated at Z.

### QUESTION 5

The following represents a generalised chromatogram from an HPLC system:



Which of the following statements is true about the above chromatogram?

- A Substance A has the highest retention time.
- B Substance C has the highest retention time.
- C Substance B has the lowest retention time.
- D Substance C has the lowest retention time.


### QUESTION 6

Suggest reasons why analysis of vitamin C (ascorbic acid) is best carried out with an HPLC system and not gas chromatography.

**Solution**

### QUESTION 7


The following is a chromatogram of the separation of five amino acids:



- Name the amino acid, which is most strongly adsorbed into the stationary phase.
- Name the amino acid with the highest affinity for the mobile phase.
- Suggest what change or changes would be necessary to alter the order and the retention times of the amino acids.

### QUESTION 8

The following is a series of chromatograms depicting the peak heights and retention times of a caffeine standard and a sample from 100ml of instant coffee. The coffee sample was prepared by dissolving a level teaspoon of instant coffee in a standard 250mL cup.



- (a) How can it be deduced that the sample peak is also caffeine?
- (b) By comparing the peak heights, estimate the amount of caffeine in the sample.
- (c) Calculate the amount of caffeine in the original 250mL cup
- (d) Suggest a more accurate way of calculating the exact amount of caffeine from the peaks.

## SOLUTIONS

**QUESTION 1** Answer is B

**QUESTION 2** Answer is A

**QUESTION 3** Answer is C

**QUESTION 4** Answer is D

**QUESTION 5** Answer is B

### QUESTION 6

Ascorbic acid is very heat sensitive and HPLC is carried out at room temperature, therefore there is less of a chance of decomposition. With GC, the sample is vaporised and also the column is heated, ascorbic acid would surely decompose at such temperature.

### QUESTION 7

- (a) Cystine
- (b) Leucine
- (c) Changes in solvent and/or column would change the retention time and even the order of exit from the column.

### QUESTION 8

- (a) The retention time is around the same. \*on a proper chromatogram it would be more evident.
- (b) If 360mg = 4 units and sample = 1 unit. Caffeine in sample = 90mg
- (c) 90mg/100mL  
 $\therefore \text{mg}/250\text{mL} = 90 \times 2.5 = 225\text{mg}$
- (d) Especially if peaks are wide, a more accurate method is to calculate the area under the peaks.